Copper(II) hydroxide carbonate

CuCO3. Cu (041) 2

Copper(II) hydroxide carbonate is a blue-green, fine crystalline powder insoluble in cold water. In hot water it decomposes into CuO, CO_2 and H_2O .

Copper(II) hydroxide carbonate is obtained according to the reaction equation:

$$2 \text{ CuSO}_4 + 2 \text{ Na}_2 \text{CO}_3 + \text{H}_2 \text{O} \longrightarrow \text{CuCO}_3 \cdot \text{Cu(OH)}_2 + \text{CO}_2 + 2 \text{ Na}_2 \text{SO}_4$$

Realization:

In $100 \, \text{cm}^3$ of water dissolve $15 \, \text{g}$ of $\text{CuSO}_4 \cdot 5 \, \text{H}_2\text{O}$ and filter into a $250 \, \text{cm}^3$ beaker. In another beaker, dissolve $18 \, \text{g}$ of $\text{Na}_2\text{CO}_3 \cdot 10 \, \text{H}_2\text{O}$. Heat both solutions to $50\text{-}60 \, ^\circ\text{C}$. Then pour the warm sodium carbonate solution into the warm copper(II) sulfate solution, stirring intensively. Check with indicator paper that the solution above the precipitate is alkaline. If not, add some more sodium carbonate solution. Heat the precipitate and mix it until $\text{CuCO}_3 \cdot \text{Cu}(\text{OH})_2$ becomes granular. Then filter and wash with warm water to remove the sulfate ions. After drying the obtained compound at $50 \, ^\circ\text{C}$, weigh it and calculate the process efficiency.

Equipment:

- · Graduated cylinder for 250 cm³,
- Beakers: 500 cm³, 2x250 cm³,
- · Heating set,
- · Filtration set,
- Stirring rod,
- · Test tubes.

Reagents:

- Copper(II) sulfate pentahydrate CuSO₄ 5 H₂O 15 g,
- · Sodium carbonate decahydrate Na₂CO₃ 10 H₂O 18 g,
- · Chloric acid 1M,
- · Ammonia solution 10%,
- · Phenolphthalein or universal indicator papers.

