

Buffer solutions

1. Experimental procedures

Compounds of group A	Compounds of group B				
0.1 M NaCl pH = 7.00	$0,1 \text{ M NH}_4 \text{ Cl}$ i $0.1 \text{ M NH}_3 \text{ aq}.$	pK _a = 9.2			
$0.1 \text{ M Na}_2 \text{ SO}_4 \text{ pH} = 7.65$	0.1 M CH ₃ COOH i 0.1 M CH ₃ COONa	$pK_a = 4.8$			
	0.1 M Na ₂ HPO ₄ i 0.1 M KH ₂ PO ₄	$pK_a = 7.2$			

1.1. The effect of addition of acid or base on the pH value of buffer and non-buffer solution

Take 30 cm³ of the solution A (depending on the choice) into the beaker in which the magnetic dipole is located and measure the pH. Then start the magnetic stirrer, dip the electrode and add 1 cm³ of 0.1 M NaOH solution (or 0.1 M HCl depending on the choice) from the burette - simultaneously taking the pH meter readings after each portion of the reagent added. The same should be done with a 30 cm³ sample, according to the choice of solution B (components mixed in a 1: 1 ratio). Write the results in Table 2.1.

1.2. Influence of mutual molar ratio of components of buffer mixture on its pH value

Use a volumetric pipet to dispense respectively 5, 10, 15, 20 cm³ of the first component of compounds of group B into a 25 cm³ flask, and make up to the mark with a second component of compounds of group B. After thorough mixing of the contents of the flask, measure the pH value. Write the results in Table 2.3.





Fundusze Europejskie Wiedza Edukacja Rozwój

Rzeczpospolita Polska Unia Europejska Europejski Fundusz Społeczny



2 Development of results

Table 2.1.

Volume of acid or base added [cm ³]				
pH salt solution (A)				
pH buffer solution (B)				

Table 2.2.

Concentration the I component of the buffer mixture [mol/dm ³]			
Concentration the II component of the buffer mixture [mol/dm ³]			
pH buffer solution			

- present results from experiment 1.1. in the form of a graph showing the dependence of the pH value on the amount of acid added (a titration curve).
- read from the titration curve how much acid (base) should be added to individual solutions to change the pH by one unit
- calculate what concentrations of strong acid (base) will correspond to pH values equal to the pH of buffer solutions (B).
- calculate the buffer capacity of the test solution from group B.

3. Conclusions

- 4. The scope of material
 - buffer capacity,
 - buffer solutions,
 - preparation of buffer solutions.

5. References

- M. D. Joesten, J. L. Wood, *World of Chemistry*, second edition, Thomson, USA 1996
- D. W. Oxtoby, N. H. Nachtrieb, *Principles of modern Chemistry*, Saunders College Publishing, USA 1996









